Angela Antoni
INDUSTRIE SpA

M O R E T H A N T E C H N O L O G Y

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advanced cold technology in the research®

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Angelantoni Industrie... realized the most complete range of blood bank refrigerators for small, medium and large Blood Centers, Hospitals, Nursing Homes, laboratories operating-theater, etc. The models shown in this leaflet have capacities from 36 to 1540 bags. Particular attention has been paid to the operating and control system of the blood banks in order to fulfill the present standards, installing a centralized electronic system in one module only, named "E 2003". This system, which is our exclusive project and manufacture, is the result of many years of experience in this sector (more than 600 blood banks sold in the last ten years), utilizing the most modern technological systems. Operation and control functions are carried out by the "E 2003" system by means of a microprocessing intelligent circuit.

Control functions are:

- * blood bank internal temperature check
- * evaporator temperature check in the defrosting phase
- * control panel push button check

On the basis of the set data, the system controls the following:

- * requested temperature regulation
- * over/under temperature alarm system regulation, independent from the control system
- * blood bank on/off
- * regulation of the automatic defrosting system.

The operating functions are electrically carried out by the power relays controlling various electric parts. Furthermore, the system is fitted for the repetition of the remote alarms which can be connected to the most frequented rooms, according to the customer needs, such as: emergency wards, reception, porter's lodge, etc. The system works with its own 6V accumulator automatic recharge, ensuring every alarm indication in case of a power failure for long.

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FRAME

zinc plated steel self bearing body, prepainted in white for "MINI" and "BB" models; "BBP" models have zinc plated prepainted white steel prefabricated modular panels. Light blue top band. AISI 304 stainless steel internal coating for monolithic BB models.

THERMAL INSULATION

obtained by injected polyurethane foam under pressure ("sandwich" technique).

SERVICE DOORS

realized with special anodized aluminium profile with triple under vacuum glass and anti-condensing electric heaters to avoid frost formation.

Magnetic lock system with odorless rubber gaskets.

Door handles in chromate alloy with plastic profiles.

Lock with external key.

DRAWERS

to contain the bags, they are realized in stainless steel with perforated surface for a proper internal air circulation. The drawers are complete with stainless steel sliding guides, having plastic wheels. The frontscreen of the drawers is in plexiglass in order to allow immediate view of the content of the drawers and to avoid warm air coming inside while the door is opened. Every drawer has a proper partition screen in transparent PVC to hold the bags in an inclined position for easy label reading. The partition screens can be adjusted or removed for different use. On the frontscreen of each drawer special plastified label holders are installed to indicate the drawer content.

COOLING SYSTEM

completely sealed located in the lower part for the "MINI" models and located in the upper part for the other models, with a particularly noiseless hermetic compressor and air cooled condenser. Copper tube evaporator with aluminium fins and helicoidal electric fans.

INTERNAL LIGHTING

by waterproof glow lamp or fluorescent tubes with switch positioned on the control panel.

STANDARD FEATURES

- * internal ventilation
- * adjustable feet for monolithic models
- * water condensing collecting tank, with special electric heaters for water evaporation in monolithic models.
- * "E 2002" model microprocessor with electronic control and digital display.
- * over/under temperature alarm system with automatic recharge accumulator.
- * predisposition for the repetition of the remote alarm lamps.
- * electric defrosting system with electronic control.
- * battery electric temperature chart recorder with aluminium frame, door with control window and key lock, seven day cycles and recording pen point with regulation control.
- * no. 55 spare diagrams (for one year running) and no. 2 spare pen points.
- * protection fuses and other accessories according to CEI 66/5 standards.

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"E2002" microprocessor electronic
operating and control system
with digital display

divisori interni in PVC
per sacche (optional)
PVC partition screens
for bags

cassetti in acciaio inox
sliding drawers in stainless steel

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aluminium profile door
with triple glass
anti-condensing system

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




The blood bank refrigerators of the "MINI" series can be executed in a special version for installation in car blood banks, blood collecting vehicles, special blood or organ transport vans etc. The feature of this equipment is the special power supply which can be 12 V D.C. or double 220V A.C. In this last solution it is possible to use the blood bank refrigerator both while the van is parked (using A.C.) and during transport (using D.C.). A special fixing system allows the installation of the blood bank refrigerator to the van avoiding above all the vibrations during the transport. Our technical and trade departments are at your complete disposal for further information.

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
The picture has been granted by CARROZZERIA AUTODROMO, MODENA

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



Angelantoni Industrie S.p.A. in Massa Martana (Perugia, Italy) extend over an area of 80,000 square meters (more than 15,000 square meters are covered by factories and offices). Massa Martana is situated in Umbria, a region rich in art, history and traditions. No location could be more appropriate, **Angelantoni** look back to the past to better understand and anticipate the future.

That is why **Angelantoni** are the most complete and diversified European group for advanced cold technology in the research.

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biomedical division also produces:

- ❑ **H E M**  **S A F E** PC-controlled blood bank refrigerator
- ❑ -40°C and -85°C Polar freezers
- ❑ reach-in refrigerators
- ❑ -135°C / -152°C freezers (with compressors) and -180°C freezers (with LN₂ system)
- ❑ plasma shock freezers
- ❑ prefabricated cold and hot chambers
- ❑ automatic ice flakers and ice cube makers
- ❑ mortuary refrigerators and autopsy tables
- ❑ chambers for pharmaceutical stability tests and plant growth chambers (sun light)

Angelantoni
INDUSTRIE S.p.A.

laboratory and biomedical division

loc. cimacolle, 464 - 06056 massa martana (pg) - italy
tel. (++39)075.8955.1 (a.r.) - fax (++39)075.8955.312

e-mail: biomedical@angelantoni.it

<http://www.angelantoni.it>

M O R E T H A N T E C H N O L O G Y