



Angela Antoni

MCT 700

THERMOSTATIC CHAMBER



laboratory and biomedical division

TECHNICAL DATA

External dimensions (mm)	720 x 880 x 1995 (LxDxH)
Internal dimensions (mm)	600 x 735 x 1485 (LxDxH)
Internal volume (l)	650
Operative temperature (°C)	+10/ +50
Accuracy - Space	±1
Accuracy - Time	±0,5
Refrigerant gas	R134a
Electrical input	230V / 50Hz / 1 + G
Absorbed current (A)	6
Max. used power (kW)	1,4
Noise level (dB(A))	<65
Thermal dissipation (kcal/h)	400
Overall weight with standard packaging (kg)	250

